


# PD-L1 Expression on Circulating Antigen Presenting Macrophages in Blood Predicts PFS & OS in an Array of Metastatic Cancer Types Treated with PD-L1/PD-1 Immunotherapies

Daniel L. Adams<sup>1</sup>, Pablo Lopez<sup>2</sup>, Steven H. Lin<sup>2</sup>, Thai Ho<sup>3</sup>, Kirby P Gardner<sup>1,4</sup>, Cha-Mei Tang<sup>5</sup>

<sup>1</sup>Creatv MicroTech, Monmouth Junction, NJ 08852, <sup>2</sup>MD Anderson Cancer Center, Houston, TX 77030, <sup>3</sup>Division of Hematology and Medical Oncology, Mayo Clinic Cancer Center, Scottsdale, AZ 85259, <sup>4</sup>Rutgers University, School of Graduate Studies, New Brunswick, NJ 08901, <sup>5</sup>Creatv MicroTech, Rockville, MD 20850

## ABSTRACT

Cancer Associated Macrophage-Like cells (CAMLs) are a specific type of antigen presenting circulating stromal cells found in the blood of cancer patients (**Figure 1**). Recently, PD-L1 expression on CAMLs was described in a variety of cancers, appearing to parallel the inflammatory PD-L1 state of the tumor microenvironment, and acting as a possible predictive biomarker for PD-L1/PD-1 immunotherapeutics (IMTs). To determine if the predictive value of CAMLs, we monitored CAML PD-L1 from n=54 patients in an array of solid tumor subtypes, prior to induction of a new line of systemic PD-L1/PD-1 immunotherapy, to evaluate CAML's ability to predict response by progression free survival (PFS) and overall survival (OS).




**Figure 1.** Example of CTC isolated with a CAML. CTCs are Cytokeratin pos. (green) & CD45/CD14 neg. CAMLs are CD45/CD14 pos. (purple) & maybe pos. for Cytokeratin. White blood cells are smaller CD45/CD14 pos.


**Table 1. Demographic Table**

Median Age (Range)	66 (34-81)
<b>Gender</b>	
Male	23 (43%)
Female	31 (57%)
<b>Race</b>	
Caucasian	47 (87%)
AA	3 (6%)
Other	3 (6%)
<b>ECOG</b>	
0	39 (72%)
≥1	15 (28%)
<b>Node positive</b>	
0	16 (30%)
≥1	33 (61%)
unknown	5 (9%)
<b>Number Prior Lines of Therapies</b>	
1	24 (44%)
2	11 (20%)
≥3	19 (35%)
<b>Cancer Types</b>	
Lung	24 (44%)
Breast	20 (37%)
RCC	10 (19%)
<b>Immunotherapy Type</b>	
Pembrolizumab	35 (65%)
Nivolumab	14 (26%)
Atezolizumab	5 (9%)
<b>Number of Metastatic Sites</b>	
1	36 (67%)
≥2	18 (33%)

**Figure 2. PFS of CAML PD-L1**



**Figure 3. OS of CAML PD-L1**



## MATERIALS & METHODS

Three single blind multi-year prospective studies were run to evaluate PD-L1 CAML expression prior to induction of PD-L1 or PD-1 based immunotherapy in metastatic Lung Cancer (LC) (n=24), metastatic Breast Cancer (BC) (n=20) or metastatic Renal Cell Carcinoma (RCC) (n=10) to evaluate patient's PFS and OS (**Table 1**). We recruited patients with pathologically confirmed disease, progressing by PET/CT after at least one failed prior line of therapy (**Table 1**). Blood samples (7.5 mL) were taken prior to start of new PD-L1/PD-1 therapeutic regimes, i.e. pembrolizumab (n=35), or nivolumab (n=14), or atezolizumab (n=5). Whole blood was filtered by CellSieve microfilters to collect CAMLs. CAMLs were identified by CD45/CD14 positivity and stained for PD-L1. CAML's PD-L1 was scored as a binary, high or low/negative. PFS & OS hazard ratios (HRs) were evaluated by censored univariate and multivariate analysis at 24 months (**Table 2**).

## RESULTS

- CAMLs were found in 98% of all tested samples.
- Pts with high CAML PD-L1 expression treated with IMT had significantly better PFS (HR=6.7, p<0.0001) independent of other clinical parameters (**Figure 2**)
- Pts with high CAML PD-L1 expression treated with IMT had significantly better OS (HR=4.1, p=0.0013) independent of other clinical parameters (**Figure 3**)

## FUNDING SOURCE

This work was supported by a NIH grant R43CA206840. The U.S. Army Research Office (ARO) and Defense Advanced Research Projects Agency (DARPA) (W911NF-14-C-0098). The content of the information does not necessarily reflect the position or the policy of the US Government.

**Table 2. Univariate and Multivariate Table**

Variable	Uni PFS (p value)	Multi PFS (p value)	Uni OS (p value)	Multi OS (p value)
<b>CAML PD-L1 high vs Low</b>	<b>&lt;0.0001</b>	<b>&lt;0.0001</b>	<b>0.0013</b>	<b>0.0025</b>
<b>Race</b>	0.2005		0.6736	
<b>Male vs Female</b>	0.4395		0.1464	
<b>Age &lt;65 vs ≥65</b>	0.4966		0.0128	0.0158
<b>ECOG 0 vs ≥1</b>	0.6089		0.6861	
<b>Cancer Type</b>	0.7356		0.9013	
<b>Node Negative vs Positive</b>	0.8945		0.7565	
<b># Prior Lines of Therapy 1 vs ≥2</b>	0.7098		0.9895	
<b># of Metastatic sites 1 vs ≥2</b>	0.2601		0.5507	
<b>Immunotherapy type</b>	0.7795		0.4566	

## CONCLUSIONS

- These data suggests that PD-L1 expression in circulating CAMLs appears to independently predict pts with increased benefit to PD-L1/PD-1 IMTs in an array of tumor types (**Table 2**)
- Larger and more refined studies are needed to validate these findings.

## REFERENCES

1. Adams DL, et al "Circulating giant macrophages as a potential biomarker of solid tumors." *PNAS*, 111(9):3514-3519. (2014)
2. Cristofanilli M, "Liquid Biopsies in Solid Tumors" *Springer Intl Publish.* (2017)
3. Adams DL, et al. "Sequential tracking of PD-L1 expression and RAD50 induction in circulating tumor and stromal cells of lung cancer patients undergoing radiotherapy" *Clin Can Res*, 23(19): 5948-5958. (2017)